

## UREA/BUN – UV

Diagnostic reagent for the in-vitro quantitative determination of urea in human serum, plasma and urine on both manual and automated systems.


|                         |          |                         |          |
|-------------------------|----------|-------------------------|----------|
| REF: BS.1/UU02.050.0100 | 100 test | REF: BS.1/UU02.125.0250 | 250 test |
| REF: BS.1/UU02.100.0200 | 200 test |                         |          |

### CLINICAL SIGNIFICANCE

Urea is the major end product of protein nitrogen metabolism. It is synthesized by a series of reactions in the liver called the urea cycle. In the urea cycle, ammonia is converted to urea, which is carried by blood to the kidneys for elimination from the body. The circulating levels of urea depend upon protein intake, protein catabolism and kidney function. Elevated urea levels can occur due to renal impairment or in some diseases such as congestive heart failure, diabetes, infection, or during different liver diseases. Determination of blood urea nitrogen (BUN) is the most widely used screening test for renal function together with serum creatinine. Serum creatinine is another metabolic waste product freely filtered by the glomerulus, but does not undergo tubular reabsorption. Its steady rate of elimination is frequently used to generate an index or ratio with BUN values for normalized evaluations.

### METHOD PRINCIPLE

Urea in the sample is hydrolyzed enzymatically into ammonia (NH<sub>3</sub>) and carbon dioxide (CO<sub>2</sub>). Ammonia ions formed react with α-ketoglutarate in a reaction catalysed by glutamate dehydrogenase (GLDH) with simultaneous oxidation of NADH to NAD<sup>+</sup>:



The rate decrease in the NADH concentration is directly proportional to the urea concentration in the specimen. It is determined by measuring the absorbance at 340 nm.

### REAGENT COMPOSITION

|   |  |
|---|--|
| <b>R1: Standard</b><br>UREA<br>BUN  | 50 mg/dl<br>23.4 mg/dl                                 |
| <b>R2: BUFFER REAGENT</b><br>Tris Buffer ( pH 8.5)<br>α-Ketoglutarate<br>GLDH<br>Urease<br>Sodium azide | 50 mmol/L<br>10 mmol/L<br>8 KU/L<br>5 KU/L<br>8 mmol/L |
| <b>R3: STARTER REAGENT</b><br>Sodium azide<br>NADH  | > 0.20 mmol/L<br>8 mmol/L                              |

### PRECAUTIONS AND WARNINGS

Reagent to be handled by entitled and professionally educated person.  
Good Laboratories practices using appropriate precautions should be followed in:

- Wearing personnel protective equipment (overall, gloves, glasses,).
- Do not pipette by mouth.

- In case of contact with eyes or skin; rinse immediately with plenty of soap and water. In case of severe injuries; seek medical advice immediately.
- Respect country requirement for waste disposal.  
**S56:** dispose of this material and its container at hazardous or special waste collection point.  
**S57:** use appropriate container to avoid environmental contamination.  
**S61:** avoid release in environment.

For further information, refer to the UREA/BUN reagent material safety data sheet.

### REAGENT PREPARATION, STORAGE AND STABILITY

**BioScien** working solution is prepared by mixing 4 volumes of R2: Buffer and 1 volume of R3: Starter eg. 400 µl R2 +100 µl R3. Urea/BUN reagents are supplied ready-to-use and stable up to the expiry date labeled on the bottles when properly stored refrigerated at 2–8°C.

Working solution is stable for 1 month at 2 – 8 °C  
or 8 days at 15 – 25 °C.

#### Deterioration

The **BioScien** Urea/BUN reagent is normally clear, reagent turbidity or control values out of the assigned range may be an indication of reagent deterioration.

### SPECIMEN COLLECTION AND PRESERVATION

#### Serum or plasma

No special preparation of the patient is required. Ensure non haemolyzed serum or plasma are used. The only acceptable anticoagulants are heparin, EDTA and fluoride, avoid ammonium which interfere with the assay.

Stability: 1 days at 15 -25°C; 7 days at 2 -8°C; 1 month frozen at -25°C.

#### Urine

Urine samples are prediluted 1: 50 with ammonium free water prior to assay.

Stability: 1 days at 15 -25°C; 7 days at 2 -8°C; 1 month frozen at -25°C.

### SYSTEM PARAMETERS

|   |  |
|---|--|
| Wavelength  | 340 nm   |
| Optical path  | 1 cm   |
| Assay type  | Fixed Rate   |
| Direction   | Decrease   |
| Sample Reagent Ratio<br>e.g.: Reagent volume<br>Sample volume               | 1:100<br>1 ml<br>10 µl                               |
| <b>Test reading time</b><br>First read time<br>delay time<br>last read time | 90 seconds<br>30 seconds<br>60 seconds<br>90 seconds |
| Temperature   | 37°C   |
| Incubation time   | zero   |
| Zero adjustment   | Against air  |

|                      |                             |
|----------------------|-----------------------------|
| Reagent Blank Limits | Low 1.00 AU<br>High 2.00 AU |
| Sensitivity          | 0.9 mg/dL (0.15 mmol/L)     |
| Linearity            | 200 mg/dL (33.2 mmol/L)     |

## EQUIPMENT REQUIRED NOT PROVIDED

- Sterile Syringe
- Analytical tubes, automatic pipet
- Centrifuge and spectrophotometer

## ASSAY PROCEDURE

|                  | Standard | Specimen |
|------------------|----------|----------|
| Working Solution | 1.0 ml   | 1.0 ml   |
| Standard         | 10 µl    |          |
| Specimen         |          | 10 µl    |

Mix, and after 30 seconds read the absorbance A1 of the standard or specimen. Exactly 1 minute later, read the absorbance A2 of standard or specimen.

## CALCULATION

### Serum or Plasma:

concentration (mg/dl) =  $\frac{(A2-A1) \text{ specimen} \times n}{(A2-A1) \text{ standard}}$

where n = 107 mg/dl

Urine urea concentration is determined by multiplying the result by the dilution factor (50).

## QUALITY CONTROL

To ensure adequate quality control, it is recommended that normal and abnormal commercial control serum of known concentrations included in each run. If control values are found outside the defined range, check the instrument calibration, and reagent for problems. If control still out of range please contact **BioScien** technical support.

## PERFORMANCE CHARACTERISTICS

| Precision  | Within run<br>(Repeatability) |            | Run to run<br>(Reproducibility) |            |
|------------|-------------------------------|------------|---------------------------------|------------|
|            | Normal level                  | High level | Normal level                    | High level |
| n          | 20                            | 20         | 20                              | 20         |
| Mean mg/dl | 45                            | 150        | 47                              | 153        |
| SD.        | 0.7                           | 2.7        | 0.82                            | 2.81       |
| CV. %      | 1.5                           | 1.95       | 1.63                            | 2.15       |

The results of the performance characteristics depend on the analyzer used.

### Accuracy (Methods Comparison)

Result obtained from **BioScien** Urea/BUN reagent compared with commercial reagent of the same methodology performed on 20 human sera give a correlation of 0.9.

### Sensitivity

When run as recommended, the minimum detection limit of the assay is 0.9 mg/dL (0.15 mmol/L).

### Linearity

The reaction is linear up to concentration of 200 mg/dl (33.2mmol/L). Specimens showing higher concentration should be diluted 1+2 using physiological saline and repeat the assay (resultx3).

## INTERFERING SUBSTANCES

### Haemolysis

No significant interference from Erythrocyte contamination.

### Icterus

No significant interference.

### lipemia

Lipemic specimens interfere with the method of Berthelot.

### Others

Ammonium heparin should not be used as anticoagulants. Ammonium ions should be avoided since it may cause erroneously elevated results. Color development in the Berthelot reaction is suppressed by amines, thiols, steroids and ascorbic acid.

### Reducing Substances

Color development in the Berthelot reaction is suppressed by amines, thiols, steroids and ascorbic acid.

## EXPECTED VALUES

| Serum and plasma  | mg/dl          | [mmol/L]            |
|-------------------|----------------|---------------------|
| <b>Urea:</b>      |                |                     |
| Children          | 11-39          | [1.8-6.4]           |
| Adults < 65 years | 15 -50         | [2.5-8.33]          |
| Adults > 65 years | ≤ 70           | [≤ 11.66]           |
| <b>BUN</b>        |                |                     |
| Children          | 5-18           | [0.84-2.99]         |
| Adults < 65 years | 7-23.5         | [1.16-3.89]         |
| Adults > 65 years | 7-32.9         | [< 5.44]            |
| <b>Urine</b>      | <b>g/24hrs</b> | <b>[mmol/24hrs]</b> |
| <b>Urea</b>       | 20-35          | 20-35               |
| <b>BUN</b>        | [330-580]      | [330-580]           |

## DYNAMIC RANGE

0.9 - 200 mg/dl (0.15 – 33.2 mmol/L).

## REFERENCES

1. Tietz N. W., textbook of clinical chemistry. Burtis CA, Ashwood ER, Saunders W.B. 3rd Edition, 1999 p 1239-1241
2. Patton, C. J., Crouch, S. R., Anal. Chem., 49, 464-469 (1977)
3. Shephard MD, Mezzachi RD: Clin Biochem Revs, 4:61-7, 1983.
4. Laboratory reference values. Urea nitrogen (BUN). Rochester, Minn.: Mayo Foundation for Medical Education and Research; Nov. 2010.
5. Tiffany to, jansen JM, Burtis CA,Overton JB, Scott CD. Enzymatic Kinetic Rate and end Point analyses of Substrate, By USE of A Gemsac fast analyzer. ClinChem. 1972; 18:829-840.

| SYMBOLS IN PRODUCT LABELLING |                             |                                  |
|------------------------------|-----------------------------|----------------------------------|
| <b>IVD</b>                   | For in-vitro diagnostic use | Number of <n> test in the pack   |
| <b>LOT</b>                   | Batch Code/Lot number       | Caution                          |
| <b>REF</b>                   | Catalogue Number            | Do not use if package is damaged |
|                              | Temperature Limitation      | Consult Instruction for use      |
|                              | Expiration Date             |                                  |
|                              | Manufactured by             |                                  |



**Medical Device Safety Service**  
**MDSS GmbH**  
 Schiffgraben 41  
 30175 Hannover, Germany

