

BioScien PT (ISI :1.00) Concentration and INR Table

Seconds										%	INR
Control	11	11.5	12	12.5	13	13.5	14	14.5	15.00	100	1.00
	11.55	12.08	12.60	13.13	13.65	14.18	14.70	15.23	15.75	95.00	1.05
	11.88	12.42	12.96	13.50	14.04	14.58	15.12	15.66	16.20	92.50	1.08
	12.10	12.65	13.20	13.75	14.30	14.85	15.40	15.95	16.50	90.00	1.10
	12.76	13.34	13.92	14.50	15.08	15.66	16.24	16.82	17.40	85.00	1.16
	13.09	13.69	14.28	14.88	15.47	16.07	16.66	17.26	17.85	82.50	1.19
	13.42	14.03	14.64	15.25	15.86	16.47	17.08	17.69	18.30	80.00	1.22
	13.75	14.38	15.00	15.63	16.25	16.88	17.50	18.13	18.75	77.50	1.25
	14.19	14.84	15.48	16.13	16.77	17.42	18.06	18.71	19.35	75.00	1.29
	15.07	15.76	16.44	17.13	17.81	18.50	19.18	19.87	20.55	70.00	1.37
	15.51	16.22	16.92	17.63	18.33	19.04	19.74	20.45	21.15	67.50	1.41
	16.06	16.79	17.52	18.25	18.98	19.71	20.44	21.17	21.90	65.00	1.46
	16.50	17.25	18.00	18.75	19.50	20.25	21.00	21.75	22.50	62.50	1.50
	17.27	18.06	18.84	19.63	20.41	21.20	21.98	22.77	23.55	60.00	1.57
	17.93	18.75	19.56	20.38	21.19	22.01	22.82	23.64	24.45	57.50	1.63
	18.59	19.44	20.28	21.13	21.97	22.82	23.66	24.51	25.35	55.00	1.69
	19.25	20.13	21.00	21.88	22.75	23.63	24.50	25.38	26.25	52.50	1.75
	20.13	21.05	21.96	22.88	23.79	24.71	25.62	26.54	27.45	50.00	1.83
	21.01	21.97	22.92	23.88	24.83	25.79	26.74	27.70	28.65	47.50	1.91
	22.00	23.00	24.00	25.00	26.00	27.00	28.00	29.00	30.00	45.00	2.00
	23.10	24.15	25.20	26.25	27.30	28.35	29.40	30.45	31.50	42.50	2.10
	24.31	25.42	26.52	27.63	28.73	29.84	30.94	32.05	33.15	40.00	2.21
	27.17	28.41	29.64	30.88	32.11	33.35	34.58	35.82	37.05	35.00	2.47
	29.04	30.36	31.68	33.00	34.32	35.64	36.96	38.28	39.60	32.50	2.64
	31.13	32.55	33.96	35.38	36.79	38.21	39.62	41.04	42.45	30.00	2.83
	33.55	35.08	36.60	38.13	39.65	41.18	42.70	44.23	45.75	27.00	3.05
	36.30	37.95	39.60	41.25	42.90	44.55	46.20	47.85	49.50	25.00	3.30
	39.82	41.63	43.44	45.25	47.06	48.87	50.68	52.49	54.30	22.50	3.62
	44.00	46.00	48.00	50.00	52.00	54.00	56.00	58.00	60.00	20.00	4.00
	56.21	58.77	61.32	63.88	66.43	68.99	71.54	74.10	76.65	15.00	5.11
	65.78	68.77	71.76	74.75	77.74	80.73	83.72	86.71	89.70	12.50	5.98
	79.53	83.15	86.76	90.38	93.99	97.61	101.22	104.84	108.45	10.00	7.23
	101.75	106.38	111.00	115.63	120.25	124.88	129.50	134.13	138.75	7.50	9.25
	143.88	150.42	156.96	163.50	170.04	176.58	183.12	189.66	196.20	5.00	13.08

According to the WHO regulations, MNPT should be established for each lot of PT Reagents by each lab.

*MNPT: Mean Normal Plasma Time.